

Mini Surgery Diode Laser System



Wuhan Gigaa Optronics Technology Co.,Ltd.

www.gigaalaser.com

Mini Surgery Diode Laser System

SOFT TISSUE SURGERY LASER

The new Gbox portable surgery diode laser system from GIGAA is the best value on the market today. Small size body case with 15w power, this mini surgery laser system has a wavelength of 810nm & 980nm for choice; it can be used in the continuous as well as pulsed mode with a contact or non contact handpiece.

Features

- New friendly used protocol
- Big color touch screen
- Output power calibration at the end of fibers
- Preprogrammed and customizable parameters
- Compact & lightweight, easy to transport
- Quick connector battery box
- Highly integrated and modular design
- New design inside to make maintenance easy and low cost
- High quality and reliability



Main Clinical Applications


- Dentistry
- Lipolysis (Laser fat burning)
- EVLT(Endovenous laser treatment)
- Vascular
- ENT
- Dermatology
- Orthopaedics
- PLDD(Percutaneous Laser Dsic Decompression)
- Neuro surgery
- Gynaecology
- Urology
- General surgery

Laser Type	GaAAs diode Laser
Model	Gbox 15A/B
Wavelength	810nm/980nm
Maximum Power	15W
Operation Mode	CW, Single or Repeat Pulse
Pulse Duration	10 μ s-3s
Repetition Rate	0.2Hz-20KHz
Fibers	Fibers Of 200 μ m ,400 μ m And 600 μ m With SMA905 Connector
Pilot Beam	Red Diode Laser of 635nm, Power<5mW
Control Mode	True Color Touch Screen
Voltage/Current Rating	110/220 VAC, 5A, 50/60 Hz
Dimensions	215(W)*245(L)*315(H) mm
Weight	4Kg
Safety Compliance	CE 0482

Accessories



Dental Handpiece



ENT Handpiece



PM Handpiece



Bleaching Handpiece



Foot Switch



Lipolysis Handpiece



Focus Handpiece



II-4-502,International Enterprise Center,
Guanggu Road,Wuhan 430074,PRC

Ph: 86-27-67848871 67848872
Fax: 86-27-67848873
Email: info@gigaalaser.com
www.gigaalaser.com

CE 0482
ISO13485.2003 93/42/EEC